

PDF Document: Guide to IEC 60601 Series & IEC 60601-1-2 EMC Requirements

Created for: Your Reference

Date: February 25, 2025

Subject: Overview of IEC 60601 Standards and Detailed Breakdown of Electromagnetic Compatibility (EMC) Requirements per IEC 60601-1-2

1. Introduction to the IEC 60601 Series

IEC 60601 is the foundational international standard for the safety and essential performance of medical electrical equipment. Compliance is mandatory for market access in most countries. The series is structured in a hierarchy:

- IEC 60601-1 (General Standard): Sets the core requirements for all medical electrical equipment (e.g., electrical shock protection, mechanical safety).
 - IEC 60601-1-X (Collateral Standards): Address horizontal topics applicable to many devices (e.g., EMC, Usability, Alarms).
 - IEC 60601-2-X (Particular Standards): Define specific requirements for particular device types (e.g., infusion pumps, defibrillators, MRI machines). A device must comply with its particular standard *in addition to* the general standard.
-

2. Focus on IEC 60601-1-2: EMC Collateral Standard

IEC 60601-1-2 is the collateral standard for Electromagnetic Compatibility (EMC). Its purpose is to ensure that a medical device:

1. Does not emit excessive electromagnetic interference (Emissions).
2. Is immune to a reasonable level of electromagnetic interference from other sources (Immunity).

Key Concept: Risk Management

The 4th edition integrates tightly with ISO 14971 (Risk Management). The manufacturer must:

- Define the intended electromagnetic environment (e.g., hospital, home).
 - Identify how disturbances could affect the device's Essential Performance.
 - Specify acceptable performance criteria during and after testing.
-

3. EMI Emissions Limits (The Device as a Source)

The device must not emit interference beyond these typical limits.

Test	Frequency Range	Limit / Range (Typical)	Standard Reference
Conducted Emissions (AC Power Lines)	150 kHz - 30 MHz	QP: 66 dB μ V to 56 dB μ V Avg: 56 dB μ V to 46 dB μ V	CISPR 11, Class B
Radiated Emissions	30 MHz - 1 GHz	QP: 40 dB μ V/m (measured at 10m)	CISPR 11, Class B

Radiated Emissions (High Freq.)	1 GHz - 6 GHz	Peak: 54 dB μ V/m Average: 44 dB μ V/m (measured at 3m)	CISPR 11
------------------------------------	---------------------	---	----------

QP = Quasi-Peak, Avg = Average

4. EMI Immunity Levels (The Device as a Victim)

The device must function correctly when subjected to these disturbances. Performance is judged against Criteria A, B, or C.

Test	Basic Immunity Level / Test Value	Key Parameters
Electrostatic Discharge (ESD)	± 8 kV (contact) ± 15 kV (air)	150pF, 330 Ω
Radiated RF Immunity	3 V/m	80 MHz - 2.7 GHz, 80% AM, 1 kHz

Electrical Fast Transient/Burst	±2 kV (power lines) ±1 kV (I/O lines)	5/50 ns, 5 kHz
Surge	±1 kV (Line-Line) ±2 kV (Line-Earth)	1.2/50 µs (voltage)
Conducted RF Immunity	3 Vrms	150 kHz - 80 MHz, 80% AM, 1 kHz
Voltage Dips & Interruptions	0%, 40%, 70%, 95% reduction	Various durations (0.5 cycle to 300 cycles)
Power Freq. Magnetic Field	3 A/m	50/60 Hz

5. Performance Criteria (Acceptance "Value")

The standard defines the required device behavior during and after immunity tests based on risk analysis:

- Criterion A: Normal operation within specification. No degradation.

- Criterion B: Temporary degradation or loss of function is allowed, but self-recovery to normal operation is required. No operator intervention needed.
- Criterion C: Temporary loss of function is allowed, requiring operator intervention (e.g., power cycle). Must be justified by risk management.

The device must maintain its Essential Performance according to the criterion specified for each test.

6. Examples of Particular Standards (IEC 60601-2-X)

A small sample of device-specific standards:

- -2-2: High frequency surgical equipment
 - -2-4: Cardiac defibrillators
 - -2-10: Nerve and muscle stimulators
 - -2-16: Haemodialysis equipment
 - -2-19: Infant incubators
 - -2-24: Infusion pumps
 - -2-25: Electrocardiographs
 - -2-37: Ultrasonic diagnostic equipment
 - -2-49: Multiparameter patient monitors
 - -2-72: Home healthcare ventilators
-

7. Important Notes & Disclaimer

- This is a summary document. The official IEC 60601-1-2 standard contains the complete, legally definitive requirements, test methods, and tables.
- Test levels can vary based on the declared intended use environment (professional vs. home healthcare).
- National deviations exist. Always check the national adoption in your target market (e.g., EN 60601-1-2 in Europe, ANSI/AAMI ES60601-1-2 in the USA).
- Compliance must be verified by a qualified test laboratory following the detailed protocols in the basic EMC standards (IEC 61000-4-series).
-

